***Peripheral Intravascular  
Lithotripsy (IVL) Catheter***

EN English .....	2
BG Bulgarian/Български .....	5
CS Czech/Čeština .....	8
DA Danish/Dansk .....	11
DE German/Deutsch .....	14
EL Greek/Ελληνικά .....	17
ES Spanish/Español .....	20
ET Estonian/Eesti .....	23
FI Finnish/Suomi .....	26
FR French/Français .....	29
HR Croatian/Hrvatski .....	32
HU Hungarian/Magyar .....	35
IS Icelandic/Íslenska .....	38
IT Italian/Italiano .....	41
LT Lithuanian/Lietuvių k. ....	44
LV Latvian/Latviski .....	47
NL Dutch/Nederlands .....	50
NO Norwegian/Norsk .....	53
PL Polish/Polski .....	56
PT Portuguese/Português .....	59
RO Romanian/Română .....	62
SK Slovak/Slovenčina .....	65
SL Slovenian/Slovenščina .....	68
SR Serbian/Srpski .....	71
SV Swedish/Svenska .....	74
TR Turkish/Türkçe .....	77

**Shockwave Intravascular Lithotripsy (IVL) System with the Shockwave L<sup>6</sup> Peripheral Intravascular Lithotripsy (IVL) Catheter**  
**Instructions for Use (IFU)**

**For use with the Shockwave Medical, Inc. IVL Generator and Connector Cable**

**Intended Use**

The intended use of the Shockwave Peripheral Intravascular Lithotripsy System is to dilate stenotic arteries by modifying the plaque present in the target artery using vessel wall modifying attributes of the device. The calcium modification aspect of the procedure using intravascular lithotripsy is a preconditioning of the stenosis by modifying the structure and compliance of calcified plaque prior to complete balloon dilatation of the stenosis using the same device.

**Indication for Use**

The Shockwave L<sup>6</sup> Peripheral IVL System is indicated for lithotripsy-enhanced, low-pressure balloon dilatation of lesions, including calcified lesions, in the peripheral vasculature, including the iliac, femoral, ilio-femoral, popliteal, and infra-popliteal arteries. This device is not intended for use in coronary, carotid, or cerebrovascular arteries.

**Target Population**

The Shockwave Medical Peripheral IVL System is intended to treat patients with calcified, stenotic peripheral arteries in patients over 18 years old who are candidates for percutaneous therapy.

**Contents: Shockwave L<sup>6</sup> Peripheral IVL Catheter (1)**

- The following balloon sizes (balloon diameter x balloon length) are available: 8.0x30mm, 9.0x30mm, 10.0x30mm, and 12.0x30mm. Each size offering has a 110cm catheter working length.
- Folded balloon diameters are: 0.091" (2.31mm) max. for 8.0mm to 9.0mm and 0.104" (2.64mm) max. for 10.0mm and 12.0mm
- 8.0mm to 9.0mm balloon is 7 F introducer sheath compatible; 10.0mm and 12.0mm balloon is 8 F introducer sheath compatible;
- 0.018" (0.46mm) guidewire compatible (OTW - 300cm wire)

**Devices Required but not supplied By Shockwave Medical, Inc.**


- 7 F or 8 F introducer sheath
- 0.018" (0.46mm) guidewire (300cm Length)
- 5"x96" (13cm x 244cm) Cable Sleeve
- Inflator

**How Supplied**

The IVL Catheter is supplied sterile via E-beam sterilization. The IVL Catheter is intended for single use only and is not intended for reuse or re-sterilization. Do not re-sterilize as this could damage the device and lead to patient injury. Do not reuse the device as this could result in cross-contamination that could result in patient injury. Carefully inspect all packaging for damage or defects prior to use. Do not use the device if any sign of damage or breach of the sterile barrier is observed as this could lead to malfunction of the device and/or injury to the patient. Store the IVL Catheter in a cool, dark, dry place. Storage of the device in extreme conditions may affect device performance and lead to patient injury.

**Device Description**

The IVL Catheter is a proprietary lithotripsy device delivered through the peripheral arterial system of the lower extremities to the site of an otherwise difficult to treat calcified stenosis. Energizing the lithotripsy device will generate acoustic pressure pulses within the target treatment site, disrupting calcium within the lesion and allowing subsequent dilation of a peripheral artery stenosis using low balloon pressure. The IVL Catheter is comprised of an integrated balloon with an array of integrated lithotripsy emitters for the localized delivery of acoustic pressure pulses. The system consists of the IVL Catheter, an IVL Connector Cable and an IVL Generator. The IVL Catheter is available in four (4) sizes: 8.0x30mm, 9.0x30mm, 10.0x30mm, and 12.0x30mm. The L<sup>6</sup> Peripheral IVL Catheter is compatible with a 7 or 8 F sheath and has a working length of 110cm. Refer to Figure 1 below for IVL Catheter components.

**Figure 1: Shockwave L<sup>6</sup> Peripheral IVL Catheter**

The Shockwave L<sup>6</sup> Peripheral IVL Catheter shaft contains an inflation lumen, a guidewire lumen, and the lithotripsy emitters. The inflation lumen is used for inflation and deflation of the balloon with 50/50 saline/contrast medium. The guidewire lumen enables the use of a 0.018" (0.46mm) guidewire to facilitate advancement of the catheter to and through the target stenosis. The system is designed as 'Over-the-wire' (OTW) with a 110cm shaft working length, so an exchange length (300cm) guidewire is indicated. The emitters are positioned within the balloon working length for delivery of acoustic pressure pulses. The balloon is located near the distal tip of the catheter. Two radiopaque marker bands within the balloon denote the length of the balloon to aid in positioning of the balloon during treatment. The balloon is designed to provide an expandable segment of known length and diameter at a specific pressure. The proximal hub has three ports: one for inflation/deflation of the balloon, one for guidewire lumen, and one for connection to the IVL Connector Cable.

**Required Devices for the IVL Procedure**

The IVL Catheter is to be used exclusively with the IVL Generator, IVL Connector Cable and its accessories. Refer to the Shockwave Medical, Inc. *IVL Generator and Connector Cable Operator's Manual* for preparation, operation, warnings and precautions, and maintenance of the IVL Generator and IVL Connector Cable.

**Shockwave L<sup>6</sup> Peripheral IVL Catheter Balloon Compliance Chart**

Pressure	8.0 x 30mm	9.0 x 30mm	10.0 x 30mm	12.0 x 30mm
atm - kPa	Ø (mm)*	Ø (mm)*	Ø (mm)*	Ø (mm)*
2 - 203	7.79	8.49	9.63	11.52
3 - 304	7.93	8.67	9.85	11.76
4 - 405	8.09	8.88	10.08	12.04
5 - 507	8.23	9.04	10.27	12.23
6 - 608	8.40	9.20	10.49	12.44

**Note:** 2 to 4 ATM is lithotripsy treatment balloon pressure, 4 ATM is nominal balloon pressure and post-treatment angioplasty pressure and 6 ATM is RBP (Rated Burst Pressure) of the balloon.

**Shockwave L<sup>6</sup> Peripheral IVL System Sequence Chart**

The following Shockwave L<sup>6</sup> Peripheral IVL System pulsing sequence must be followed during treatment. Do not utilize a pulsing sequence other than the one outlined in the IVL System Sequence Chart below. Insertion of any size L<sup>6</sup> IVL Catheter into the IVL Generator will automatically program it with the following treatment sequence:

Treatment Frequency	2 Hz (1 Pulse every 0.5 Seconds)
Maximum Number of Continuous Pulses (1 cycle)	30 Pulses
Minimum Pause Time	10 Seconds
Maximum Total Pulses Per Catheter	300 Pulses

In the event the user attempts to deliver more than the maximum number of continuous pulses allowed, the IVL Generator is designed to stop automatically. To resume pulsing, wait at least the minimum pause time before resuming therapy. The therapy button must be released and pressed again to resume therapy. For more information, refer to the *IVL Generator and Connector Cable Operator's Manual*.

If the maximum pulse count is reached as displayed on the generator, the catheter can no longer be used. If further therapy is needed, discard this catheter and obtain a new one.

**Caution: Do not exceed 180 pulses in the same treatment segment.**

**Contraindications for Use**

The IVL System is contraindicated for the following:

1. Unable to pass 0.018" (0.46mm) guidewire across the lesion.
2. This device is not intended for treatment of in-stent restenosis.
3. This device is not intended for use in coronary, carotid, or cerebrovascular arteries.

**Warnings**

1. This device is intended for single (one-time) use only. DO NOT re-sterilize and/or reuse.
2. Do not use a device past the expiration date on the label. Use of expired product may result in patient injury.
3. Always insert the IVL Connector Cable into a sterile sleeve prior to use.
4. Use only an appropriately sized balloon for the vessel to be treated.

5. Inflate the balloon according to the balloon compliance chart. Balloon pressure should not exceed the rated burst pressure (RBP).
6. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
7. Use the IVL Generator in accordance with recommended settings as stated in the Operator's Manual. Do not attempt to override the lifetime pulse limits per device as defined in the IVL System Sequence Chart.
8. This device should only be used by physicians who are familiar with interventional vascular procedures.
9. Physicians must read and understand these instructions prior to use of the device.
10. Do not use excessive force/torque when using this device as this could result in damage to the device components and patient injury.
11. Inspect all product components and packaging prior to use. Do not use the device if it or the packaging has been damaged or if sterility has been compromised. Damaged product could result in patient injury.
12. For preparation, operation, warnings and precautions, and maintenance of the IVL Generator and its accessories refer to the *IVL Generator and Connector Cable Operator's Manual*.

**Precautions**

1. Perform all device manipulations under adequate fluoroscopic guidance.
2. Use only the recommended balloon inflation medium.
3. Appropriate anticoagulant therapy should be administered by the physician.
4. Decision regarding use of distal protection should be made based on physician assessment of treatment lesion morphology.
5. Care should be taken not to kink the catheter. If kinking occurs, remove device and prepare a new device.
6. If an inability to inflate or maintain pressure occurs, remove the catheter and use a new device.
7. If the catheter appears not to deliver lithotripsy acoustic pressure pulses, remove and replace it with another catheter.
8. Precaution should be taken when handling device after exposure to patient, e.g. contact with blood. Used product is considered biohazardous material and should be disposed of properly as per hospital protocol.

**Expected Clinical Benefits**

The following data represented is from Shockwave M<sup>5+</sup>/5<sup>4</sup> Peripheral IVL Catheters for which the Shockwave L<sup>6</sup> Peripheral IVL Catheter is equivalent.

The clinical benefits of the Peripheral IVL System when used as intended for lithotripsy-enhanced, low-pressure balloon dilatation of calcified, stenotic peripheral arteries in patients who are candidates for percutaneous therapy include: (1) final residual stenosis <50%; (2) low incidence of flow-limiting dissections at the final angiographic timepoint; and (3) freedom from new-onset Major Adverse Events (MAE) through 30 days including: death, target limb revascularization and unplanned target limb amputation (above the ankle).

The Disrupt PAD III Study was a global prospective, multi-center, single blind, randomized clinical investigation of the Shockwave Peripheral Intravascular Lithotripsy (IVL) System. The objective of the randomized study was to assess the safety and effectiveness of IVL used in combination with a drug-coated balloon (DCB) versus standard percutaneous transluminal angioplasty (PTA) used in combination with a DCB to treat moderate and severely calcified femoropopliteal arteries. The study was designed to enroll a minimum of 334 and up to 400 evaluable subjects from up to 60 global sites. The primary effectiveness endpoint was Procedural Success defined as residual stenosis ≤30% without flow-limiting dissection (≥ grade D) prior to DCB or stenting, as assessed by angiographic core lab. The powered secondary endpoint was Primary Patency at 12 months defined as freedom from clinically-driven target lesion revascularization (TLR) and freedom from restenosis determined by duplex ultrasound (DUS) or angiogram <50% stenosis. Other secondary endpoints to be assessed at 30 days, 6, 12 and 24 months included major adverse events (MAE), clinically-driven target lesion revascularization (CD-TLR), ankle brachial index (ABI), Rutherford Category, and quality of life (QoL) measures (EQ-5D and Walking Impairment Questionnaire (WIQ)). A MAE was defined as the need for emergency surgical revascularization of target limb, unplanned target limb major amputation (above the ankle), symptomatic thrombus or distal emboli that required surgical, mechanical, or pharmacologic means to improve flow and extend the hospitalization, as well as perforations that required an intervention, including bail-out stenting. The primary analysis dataset was the intent-to-treat (ITT) population. Independent core

laboratories provided an unbiased assessment of all imaging utilized in the endpoint assessments and an independent clinical events committee (CEC) adjudicated all MAE, revascularizations and deaths.

**Table 1.** PAD III Randomized Clinical Trial (RCT) Primary Endpoint Data

Primary Endpoint <sup>1</sup>	IVL N=153	PTA N=153	P value
Procedural Success <sup>2</sup>	96/146 (65.8%)	67/133 (50.4%)	0.0065

1. The hypothesis is tested using Fisher's Exact test with a one-sided  $\alpha=0.025$ .  
2. Procedural success is defined as residual stenosis  $\leq 30\%$  without flow-limiting dissection ( $\geq$  grade D) prior to DCB or stenting by angiographic core lab. A total of 279 subjects had analyzable angiographic images (IVL 146 vs 133 PTA).

The Disrupt PAD III Observational Study (PAD III OS) was a global, prospective, multi-center, single-arm, registry of the Shockwave Peripheral Intravascular Lithotripsy (IVL) System. The objective of this study was to assess the real-world acute performance of IVL in the treatment of calcified, stenotic, peripheral arteries. The study was initially designed to enroll a maximum of 250 subjects from up to 60 global sites. The sponsor chose to increase the sample size to 1500 subjects to allow for large registry and sub-group analyses in multiple lesions with the use of the M<sup>3</sup> and/or S<sup>4</sup> catheters.

The primary effectiveness endpoint was Procedural Success defined as final residual stenosis  $\leq 30\%$  without flow-limiting dissection ( $\geq$  grade D) by angiographic core lab. An independent angiographic core laboratory provided an unbiased assessment of all imaging utilized in the endpoint assessment.

**Table 2.** PAD III OS Primary Endpoint Data

Procedural Success (Final) - Subject Level <sup>1</sup>	N=1262
Yes	69.4% (772/1112)
No	30.6% (340/1112)
Procedural Success (Final) - Lesion Level <sup>2</sup>	N=1531
Yes	73.0% (998/1367)
No	27.0% (369/1367)

1. Procedural success by subject is considered when all the treated lesions in a subject meet the Procedural Success criteria.  
2. Procedural Success is considered on lesions who had final residual stenosis  $\leq 30\%$  and final dissection grade below Grade D.  
*Note: 1367 is the total number of lesions with analyzable angiographic images evaluable to assess the primary endpoint results.*

The Disrupt PAD<sup>®</sup> Study was a prospective, multi-center, single-arm study, conducted in Australia, New Zealand and the US, designed to evaluate the safety and performance of the Shockwave M<sup>3</sup> IVL System to treat calcified peripheral arteries.

All subjects were eligible for assessment of the primary safety endpoint (Table 3). No subjects experienced an MAE, resulting in a 30-day MAE rate of 0.0% (0/37). Of 52 target lesions, 48 (92.3%) had evaluable angiographic images at the final timepoint for analysis of the primary effectiveness endpoint. Technical Success was achieved in 89.6% (43/48) of target lesions. The remaining five (5) lesions had residual stenosis  $> 30\%$ ; none had flow-limiting dissection at the final timepoint.

**Table 3.** PAD<sup>®</sup> 30-Day Endpoint Results

Measure	Statistics % (n/N)
<b>Primary Safety Endpoint</b> MAE Event	0.0% (0/37)
Need for emergency surgical revascularization of target limb	0.0% (0/37)
Unplanned target limb major amputation (above the ankle)	0.0% (0/37)
Symptomatic thrombus or distal emboli <sup>1</sup>	0.0% (0/37)
Perforations require intervention to resolve, including bail-out stent	0.0% (0/37)
<b>Primary Effectiveness Endpoint</b> Final Technical Success <sup>2</sup>	89.6% (43/48)

Measure	Statistics % (n/N)
1. Defined as clinical signs/symptoms of thrombus or distal emboli detected in the treated limb in the area of, or distal to, the treated lesion after the index procedure, or noted angiographically, that require surgical, mechanical, or pharmacological means to improve flow and extend hospitalization 2. Technical success: final residual stenosis $\leq 30\%$ without flow-limiting dissection ( $\geq$ Grade D) of the lesion as assessed by the angiographic core lab	

#### Adverse Effects

Possible adverse effects consistent with standard angioplasty include:

- Access site pain
- Allergic reaction to contrast medium, anticoagulant and/or antithrombotic therapy
- Arterial dissection
- Arterial perforation or rupture
- Arterial spasm
- Arteriovenous fistula
- Bleeding complications
- Death
- Emboli (air, tissue, thrombus, or atherosclerotic emboli)
- Emergency or non-emergency arterial bypass surgery
- Entry site complications
- Fracture of the guidewire or any component of the device that may or may not lead to device embolism, serious injury or surgical intervention
- Hematoma at the vascular access site(s)
- Hemorrhage
- Hypertension/Hypotension
- Infection/sepsis
- Ischemia
- Placement of a stent
- Pseudoaneurysm
- Renal failure
- Restenosis of the treated segment
- Shock/pulmonary edema
- Total occlusion of the peripheral artery
- Vascular complications which may prolong procedure and/or which may require surgical repair (conversion to open surgery)

Risks identified as unique to the device and its use:

- Allergic/immunologic reaction to the catheter material(s)
- Device malfunction or failure

#### Procedural Steps

**Caution:** Refer to the IVL Generator and Connector Cable Operator's Manual for preparation, operation, warnings and precautions, and maintenance of the IVL Generator and IVL Connector Cable.

#### Preparation

1. Prepare the insertion site using standard sterile technique.
2. Achieve vascular access using physician's preferred access site.
3. Place an appropriately sized introducer sheath.
  - a. For radial access, use preferred sheath of appropriate length to reach iliac bifurcation.
4. Select a balloon catheter size that is 1.1:1 based on Balloon Compliance Chart (above) and reference vessel diameter. The largest diameter balloon should be used if 1.1:1 sizing is not available (such as using a 8.0mm IVL Catheter in vessels with a reference diameter of greater than 8.0mm and up to 10.0mm).
5. Verify the product label matches the catheter selected in the previous step.
6. Inspect the sterile barrier and ensure it is intact.
7. Open the sterile barrier by peeling away the white flap from the clear pouch.
8. Carefully introduce the catheter aseptically to the sterile field
9. Prepare the balloon using standard technique. Fill a 20cc (mL) syringe with 5cc (mL) of 50/50 saline/contrast medium. Attach syringe to inflation port on catheter hub. Pull vacuum at least 3 times, releasing vacuum to allow the fluid to replace the air in the catheter.
10. Fill indeflator device with 10cc (mL) of 50/50 saline/contrast medium. Disconnect syringe and connect indeflator to inflation port of catheter hub ensuring no air is introduced to the system.
11. Flush the guidewire port with saline.
12. Remove the protection sheath from the catheter.
13. Wet the balloon with sterile saline.
14. Insert the IVL Connector Cable into a sterile sleeve or probe cover.
15. Remove the cap from the proximal end and attach the IVL Catheter's connector (see Figure 1) to the IVL Connector Cable.
16. Attach the other side of same IVL Connector Cable to the IVL Generator.

**Caution: Do not press the therapy button unless the balloon is filled with 50% saline / 50% contrast medium because this may damage the balloon.**

#### Delivering the IVL Catheter to the Treatment Site

1. Advance the 0.018" (0.46mm) guidewire across the treatment site.
2. Load the IVL Catheter over the exchange length (300cm) 0.018" (0.46mm) guidewire and through the sheath, and advance balloon to the treatment site.
3. Position the balloon at the treatment site using the marker bands to aid in positioning.

#### Treating the Site with Lithotripsy

1. Once the IVL Catheter is in place, record position using fluoroscopy.
2. If position is incorrect, adjust the IVL Catheter to the correct position.
3. Inflate balloon, not exceeding 4 atm, to ensure there is full apposition to the vessel wall.  
*NOTE: Lithotripsy should not be delivered if the balloon is inflated to  $>4$  atm as there is no increase in sonic output, and higher pressure during treatment can increase the risk that the balloon loses pressure.*
4. Deliver IVL treatment (up to 30 pulses) per the IVL System Sequence Chart by pressing the therapy button on the IVL Connector Cable.
5. Inflate balloon to nominal pressure per Balloon Compliance Chart (above) and record lesion response on fluoroscopy.
6. Deflate balloon and wait at least 10 seconds to re-establish blood flow.  
*NOTE: The IVL Generator is programmed to force a minimum pause time of 10 seconds following every 30 pulses delivered.*
7. Repeat steps 3, 4, 5, and 6 for additional treatment cycles until the lesion has been sufficiently dilated or if the catheter is re-positioned.
8. Additional treatment cycles can be performed if deemed necessary. If multiple inflations are required due to a lesion length greater than the IVL balloon length, the recommended balloon overlap is at least 1cm to prevent geographic miss. However, care must be taken not to exceed 180 pulses in the same treatment segment.
9. Perform a completion arteriogram to assess post-IVL treatment result.
10. Confirm that the balloon is fully deflated prior to removing the IVL Catheter.
11. Remove the IVL Catheter. If there is difficulty in removing the device through the hemostatic valve due to the lubricity, gently grasp the catheter with sterile gauze.
12. Inspect all components to ensure that the catheter is intact. If a device malfunction occurs or any defects are noted on inspection, flush the guide wire lumen and clean the outer surface of the catheter with saline, store the catheter in a sealed plastic bag, and contact Shockwave Medical, Inc. for further instructions.

**Caution: IVL Catheter once pulled out of the body should not be reinserted for additional inflation or lithotripsy treatments. Balloon can be damaged in the process.**

#### Patient Information

Physicians should instruct patients to seek medical attention immediately for signs and symptoms of decreased peripheral blood flow. There are no known limitations to normal daily activities. Patients should be instructed to comply with the medical regimen as prescribed by their physician.

#### Return of Devices

If any portion of the Shockwave Medical IVL System fails prior to or during a procedure, discontinue use and contact local representative and/or email [complaints@shockwavemedical.com](mailto:complaints@shockwavemedical.com).

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices): if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

**Patents:** [www.shockwavemedical.com/patents](http://www.shockwavemedical.com/patents)

Symbol	Definition
	Do not re-use
	Medical device
	Use by date
	Sterilized using irradiation; Single sterile barrier with protective packaging outside
	Sterilized using irradiation
	Caution
	Manufacturer
	Date of Manufacture
	Procedure Pack Producer
	Do not use if package is damaged; If sterile barrier is damaged, consult instructions for use
	Keep dry
	Authorized Representative in the European Community
	Keep away from heat
	Batch code
	Catalogue number
	Do not resterilize
	Crossing profile

Symbol	Definition
	Non-pyrogenic
	Consult instructions for use
	Contains 1 unit (Contents: 1)
	Recommended Guidewire
	Recommended Introducer Sheath
	Over-the-wire
	Balloon Diameter
	Balloon Working Length
	Catheter Working Length (Usable Length, UL)
	Conformité Européenne
	Patents. Refer to <a href="http://www.shockwavemedical.com/patents">www.shockwavemedical.com/patents</a>
	Peripheral Arterial Disease
	Indicates a carrier that contains Unique Device Identifier information.
	Importer
	Authorized Representative in Switzerland

Shockwave Medical, Inc.  
5403 Betsy Ross Drive  
Santa Clara, CA 95054, USA  
[www.shockwavemedical.com](http://www.shockwavemedical.com)

ICON (LR) Limited  
South County Business Park  
Leopardstown  
Dublin 18  
D18 X5R3  
Ireland

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